



PRESTAN®

PRESTAN Ultralite® Manikin 12-Pack

Prestan's most portable manikin, the **Ultralite® Manikin** is amazingly easy to set up for convenient and lightweight training-on-the-go! This durable manikin offers an affordable method with CPR Feedback training, with all the quality and realism you expect from Prestan.



Prestan's most portable manikin, the **Ultralite® Manikin** is now available in an **ULTRA-CONVENIENT 12-PACK!**

- > Wheels and telescoping handle make for effortless maneuverability
- > Only 39 pounds for easy lifting in and out of vehicles
- > Stackable packing system for portability

**A CLASSROOM
OF CPR TRAINING
MANIKINS IN ONE
PORTABLE PACK!**



Ultralite® Manikin Features:

- Uniquely fast and easy set-up
- Clamshell opening provides easy access
- Manikin is easy to clean
- Vinyl and latex free
- Compatible with AED Trainer Pads
- Anatomically correct design includes xyphoid process
- Available in medium and dark skin tones
- Realistic head tilt
- Visible chest rise with lung bag use
- Realistic chest compressions to a 2" depth
- Three (3) year warranty

Ultralite® Manikin 12-pack with Feedback

PP-ULM-1200M-MS

PP-ULM-1200M-DS

Ultralite® Manikin 12-pack without Feedback

PP-ULM-1200-MS

PP-ULM-1200-DS

Reorder consumables from your Prestan Authorized Distributor:

- #PP-ULB-50 (50-pack Ultralite Manikin face-shield/lung-bags)
- #PP-UFS-50 (50-pack Ultralite Manikin face-shields)

