



PRESTAN Professional AED Trainer

The Prestan Professional AED Trainer provides quality, durability and affordability like no other trainer in the market. Clear, confident voice instructions that provide easy-to-follow directions and preconfigured scenarios make training students easier. The Prestan Professional AED Trainer carries a three (3) year warranty.

Features:



- Easy to use panel includes a scenario, volume control, and language button
- Includes 5 scenarios for varied training opportunities



- Pads can be used at least 50-60 times during CPR/AED classes
- Pads are pre-connected

ITEM # PP-APAD-1 (One set of Adult Pads)
ITEM # PP-APAD-4 (Four sets of Adult Pads)
ITEM # PP-ACASE-1 (One Adult Pads Case)



ITEM # PP-PPAD-1 (One set of Pediatric Pads)
ITEM # PP-PPAD-4 (Four sets of Pediatric Pads)
ITEM # PP-PCASE-1 (One Pediatric Pads Case)
ITEM # PP-PPDC-1 (One set of Pediatric Pads w/Case)



- Pads have an embedded 'switch' which senses proper pad attachment to the manikin
- The silicone base adhesive on the pads, coupled with the foam structure, ensures easy application & removal
- Pads carry a one (1) year warranty



- Module contains the most up to date guidelines and language capability
- When guidelines change, only the module needs to be replaced
- Module contains 2 languages (English / Spanish)



Optional Remote Control:

- Simulate loose pad
- Initiate "Press Deeper" during CPR
- Pause, play or stop
- Volume, language and scenario control

ITEM #PP-AEDT-100-R (Single Remote)
ITEM #PP-AEDT-400-R (Four Remotes)



Patent Pending

ITEM # PP-AEDT-101 (Single AED Trainer)

ITEM # PP-AEDT-401 (Four AED Trainers)

- Includes one set of Adult Pads w/Case, batteries, instructions and a carrying case
- Compatible with any CPR manikin
- Can be used to train in both semi-automatic and fully-automatic versions
- Metronome and counts upward by 10's ending at 30 to help the student with chest compressions

ITEM # PP-AEDT-KIT-101 (Kit includes Single AED Trainer, Adult Pads w/Case, Pediatric Pads w/Case and Remote)

