




## PRESTAN® Professional Infant Manikin

### Infant Manikin features:

- Life-like features and weight allows for realistic infant CPR training
- A unique clamshell design that accommodates an easy-to-insert face-shield/lung-bag
- Realistic-feeling Latex-free skin that won't shrink or stretch, even with repeated cleanings
- A clicker mechanism that helps students learn to use the right amount of force in order to compress the chest to the correct depth – great for boosting students' confidence in their ability to administer adequate compressions



- A head that simulates the way an actual infant's head would move. Face/head tilt enhances the CPR learning experience and teaches students how to open the airway during CPR
- Full chest rise with lung bag ensures students administer the most accurate rescue breaths
- Face shield is secured so it won't slip and there's no need for adhesives
- Available individually or in convenient multipacks of four
- Located on the diaper is an optional CPR Rate Monitor that provides instant visual feedback on compression rates - students know they are at the correct compression rate when two green lights illuminate



**ITEM #PP-IM-400M-MS (with CPR Monitor)**  
**(PRESTAN Professional Infant Manikin 4-Pack)**

### Reorder Accessories:

- **#PP-ILB-50**  
(50-pack Manikin face-shield/lung-bags)
- **#PP-IFS-50**  
(50-pack Manikin face-shields)



**ITEM #PP-IM-100M-MS (with CPR Monitor)**  
**(PRESTAN Professional Infant Manikin Single)**