



***Building Confidence***  
**with PRESTAN**  
CPR Manikins and AED Trainers




# PRODUCT CATALOG

---

[WWW.PRESTAN.COM](http://WWW.PRESTAN.COM)



# Building Confidence with Feedback...

PRESTAN offers a full range of CPR Training Manikins and AED Trainers designed to meet the needs of all training environments - including **Portable**, **Professional**, and **Advanced** training scenarios. Whether an instructor is starting a new business or looking to upgrade their equipment, PRESTAN offers the most robust product offering in the market giving students the *confidence* and skills necessary to be ready for all lifesaving moments.



PRESTAN CPR Training Manikins have provided a feedback mechanism since their inception. All **PRESTAN Professional Manikins** and **PRESTAN Ultralite Manikins** come equipped with a CPR Monitor that provides both visual and audible feedback, helping students learn the proper compression rate and depth and allowing them to *build confidence* and skills needed to perform high quality CPR.

The **PRESTAN Professional Series 2000 Manikins** provides advanced feedback through the addition of a Bluetooth connection to the **PRESTAN CPR Feedback App** which monitors rate, depth, recoil, ventilation, hands-off time, and chest compression fraction.



The **PRESTAN Professional AED Trainer PLUS**, the **PRESTAN AED UltraTrainer**, and the **PRESTAN AED Raizor Trainer** provide a realistic training experience with clear, confident voice instructions and easy-to-follow CPR coaching feedback.



# Building Confidence with PRESTAN...



PRESTAN is committed to *building confidence* in CPR/AED training instructors and their students.

PRESTAN offers a full range of CPR Training Manikins and AED Trainers which are intuitive, durable, cost-effective, and provide immediate feedback.

PRESTAN also offers a full range of Training Supplies and Accessories including replacement lung bag/face shields, stand-alone face shields, training face masks & adaptors, manikin clothing, and kneeling pads.



## Find a PRESTAN Authorized Distributor

PRESTAN products are sold exclusively through a global network of PRESTAN Authorized Distributors.

Please visit the 'Find a Distributor' page on the PRESTAN website to locate a PRESTAN Authorized Distributor near you.



## Visit the PRESTAN Online Resource Center

PRESTAN provides an extensive set of educational tools and training resources that will give you and your students the *confidence* to be ready.

The PRESTAN Resource Center contains Product Instruction Sheets, Frequently Asked Questions, Troubleshooting & Product Help Guides, Instructional How-To & Product Overview Videos, and much more!



# Building Confidence with PRESTAN...



PRESTAN *builds confidence* in instructors by offering an industry-leading **3-year Limited Warranty** on all PRESTAN CPR Training Manikins, the Professional AED Trainer PLUS and the AED UltraTrainer. The PRESTAN AED Raizor Trainer is backed by a 1-year Limited Warranty.

To learn more about the PRESTAN Limited Warranty, please contact a PRESTAN Authorized Distributor.



PRESTAN Products, LLC is headquartered in the United States and operates in the *heart* of the country near Cleveland, Ohio.

Unlike any other branded or private labelled CPR Training Manikin available throughout the world, all PRESTAN CPR Training Manikins are **Made in the USA**.

The company was founded in 2004 and has grown from a small start-up to a global company with warehousing facilities located throughout the world.

To learn more about PRESTAN Products, please visit **[www.prestan.com](http://www.prestan.com)**.





# INNOVATIVE CPR & AED TRAINING SOLUTIONS



## PORTABLE

- Ultralite® Adult Manikin
- Infant Ultralite® Manikin
- AED UltraTrainer
- AED Raizor Trainer

## PROFESSIONAL

- Professional Adult/PRO+ Manikin
- Professional Female Manikin
- Professional Child Manikin
- Professional Infant Manikin
- Professional AED Trainer PLUS



## ADVANCED

- Professional Adult/PRO+ Series 2000 Manikin
- Professional Female Series 2000 Manikin
- Professional Child Series 2000 Manikin
- Professional Infant Series 2000 Manikin

## WHY CHOOSE PRESTAN®

- ✓ Industry Compliance
- ✓ Unmatched Realism
- ✓ Simple & Intuitive

- ✓ Robust Product Offering
- ✓ Durability You Can Trust
- ✓ Made in the USA



In combination with PRESTAN's traditional CPR feedback light system, the **Professional Adult/Female Series 2000 Manikins** offer a unique, dual-purpose training solution providing instructors with an option to train with either basic or advanced CPR feedback including compression rate, depth, recoil, ventilation, hands-off time, and chest compression fraction via a Bluetooth® enabled PRESTAN CPR Feedback App.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



### Professional Adult Series 2000, Single

PP-AM-2000-1-MS

PP-AM-2000-1-DS



### Professional Female Series 2000, Single

PP-AFM-2000-1-MS

PP-AFM-2000-1-DS



\*Mobile device not included

### Professional Adult Series 2000, 4-pack

PP-AM-2000-4-MS

PP-AM-2000-4-DS




### Professional Female Series 2000, 4-pack

PP-AFM-2000-4-MS

PP-AFM-2000-4-DS

The **PRESTAN CPR Feedback App** provides high-quality, real-time, advanced CPR feedback via a Bluetooth® connection with the PRESTAN Professional Adult/Female Series 2000 Manikin.



**Download Now**  
on Apple or Android devices



- Manikin shirts included with every manikin.



The **PRESTAN PRO+ Adult Series 2000 Manikin** offers Advanced Feedback including visual chest rise, rate, depth, recoil, ventilation (through mouth & nose), hands-off time, and chest compression fraction via a Bluetooth® enabled PRESTAN CPR Feedback App available on both Apple and Android devices.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



The PRO+ Adult Series 2000 Head design is enhancing realism with a functional nasal airway and the ability to use a NARCAN®, OPA, or NPA training device.



### PRO+ Adult Series 2000, Single

PP-AMP-2000-1-MS

PP-AMP-2000-1-DS



### PRO+ Adult Series 2000, 4-pack

PP-AMP-2000-4-MS

PP-AMP-2000-4-DS



### PRO+ Lung Bag

PP-VAPLB-50 (50-pack)



### PRO+ Adult Face Shield

PP-APFS-50 (50-pack)

The **PRO+ Adult Face Shield** has been designed to fully cover both oral and nasal airways; the standard PRESTAN Professional Adult Face Shield can also be used.

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin

The **PRESTAN Professional Child/Infant Series 2000 Manikins** combine PRESTAN's trusted realism with advanced CPR feedback metrics, including compression rate, depth, recoil, ventilation, hands-off time, and chest compression fraction.



**Professional Child Series 2000, Single**

PP-CM-2000-1-MS  
PP-CM-2000-1-DS



**Professional Infant Series 2000, Single**

PP-IM-2000-1-MS  
PP-IM-2000-1-DS



**Professional Child Series 2000, 4-pack**

PP-CM-2000-4-MS  
PP-CM-2000-4-DS



**Professional Infant Series 2000, 4-pack**

PP-IM-2000-4-MS  
PP-IM-2000-4-DS

The **Professional Child/Infant Series 2000 Lung Bags** use the same easy to install PRESTAN design with the addition of a breath sensor. The Series 2000 lung bags are required for ventilation feedback via the PRESTAN CPR Feedback app.



**Series 2000 Lung Bags**

PP-VCLB-50 (50-pack)  
PP-VILB-50 (50-pack)

The **PRESTAN CPR Feedback App** provides high-quality, real-time, advanced CPR feedback via a Bluetooth® connection.



**Download Now**  
on Apple or Android devices



- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin




# PROFESSIONAL

## PRESTAN Professional Adult Manikin



Realistic to the eye and the touch, the **Professional Adult Manikin** is unlike any other CPR training manikin on the market. Designed to provide each student with a real-life training experience along with immediate compression rate and depth feedback, PRESTAN manikins help instructors *build confidence* in their students when teaching high-quality CPR skills.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



### Professional Adult, Single

PP-AM-100M-MS

PP-AM-100M-DS



### Professional Adult, 4-pack

PP-AM-400M-MS

PP-AM-400M-DS



### Professional Adult Diversity Kit

PP-AM-400M-MSDS

(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin

# PROFESSIONAL

## PRESTAN Professional Female Manikin



The **Professional Female Manikin** is designed to transform CPR/AED training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration. Using the same trusted design as the Professional Adult Manikin, the Female Manikin is realistic to the eye and the touch, unlike any other CPR Training Manikin on the market.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



When used in combination with the PRESTAN Female Accessory (to represent a bra), Manikin Shirt, and AED Trainer, PRESTAN offers the most realistic training experience to help *build confidence* when learning these critical lifesaving skills.

### Professional Female, Single

PP-AFM-100M-MS

PP-AFM-100M-DS



### Professional Female, 4-pack

PP-AFM-400M-MS

PP-AFM-400M-DS



### Professional Female Diversity Kit

PP-AFM-400M-MSDS

(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin



# PROFESSIONAL

## PRESTAN Professional Female Manikin



Designed to support best practices in CPR education, this mix of 2 Adult and 2 Female Professional Manikins encourages balanced gender representation in the classroom and helps address disparities in CPR response between men and women.



**Professional Adult/Female, 4-pack**

PP-AM-401M-MS  
PP-AM-401M-DS



**Professional Adult/Female  
Diversity Kit, 4-pack**

PP-AM-401M-MSDS



- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin




# PROFESSIONAL PRESTAN PRO+ Manikin



The **PRO+ Adult Manikin** features a fully functioning oral and nasal airway along with an improved lung bag design for fast and easy installation. With the introduction of the PRO+ Adult Head design, PRESTAN is enhancing realism with a functional nasal airway and the ability to use a NARCAN®, NPA, or OPA training device.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



The PRO+ Adult Head includes functional nasal and oral airways and the ability to use a NARCAN®, OPA, or NPA training device.



## PRO+ Adult, Single

PP-AMP-100M-MS  
PP-AMP-100M-DS



## PRO+ Adult, 4-pack & Diversity Kit

PP-AMP-400M-MS  
PP-AMP-400M-DS  
PP-AMP-400M-MSDS

The **PRO+ Adult Head** is also sold separately, allowing you to upgrade any Professional Adult/Female Manikin.



## PRO+ Adult Head

PP-APHEAD-1-MS  
PP-APHEAD-1-DS  
PP-APHEAD-4-MS  
PP-APHEAD-4-DS



## PRO+ Lung Bag

PP-APLB-50 (50-pack)



## PRO+ Adult Face Shield

PP-APFS-50 (50-pack)

The **PRO+ Adult Face Shield** has been designed to fully cover both oral and nasal airways; the standard PRESTAN Professional Adult Face Shield can also be used.

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin



# PROFESSIONAL

## PRESTAN Professional Child/Infant Manikins



The **Professional Child Manikin** is designed to look like a child versus an adult with less prominent body structures and softer, more childlike facial features. It is smaller and thinner than the Professional Adult Manikin yet still realistic to the eye and touch.



### Professional Child, Single

PP-CM-100M-MS  
PP-CM-100M-DS

The **Professional Infant Manikin** life-like features and weight allows for a realistic infant CPR training experience. Sculpted body contours and realistic skin feel help prepare students for real-life rescues with a head designed to simulate the way an actual infant's head moves.



### Professional Infant, Single

PP-IM-100M-MS  
PP-IM-100M-DS



### Professional Child, 4-Pack

PP-CM-400M-MS  
PP-CM-400M-DS



### Professional Child Diversity Kit

PP-CM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)



### Professional Infant, 4-Pack

PP-IM-400M-MS  
PP-IM-400M-DS



### Professional Infant Diversity Kit

PP-IM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin

# PROFESSIONAL

## PRESTAN Specialty Packs



The **Professional Family Pack** combines two Professional Adult, one Child, and two Infant Manikins in one single collection.



### Professional Family Pack

PP-FM-500M-MS  
PP-FM-500M-DS

The **Professional Collection** combines a Professional Adult, Child, and Infant Manikin in one single collection.



### Professional Collection

PP-FM-300M-MS  
PP-FM-300M-DS

The **Professional Adult Jaw Thrust Manikin** gives instructors the option to teach with the jaw thrust or head-tilt/chin-lift maneuver allowing the student to grasp the manikin's lower jaw and lift with both hands displacing the jaw of the manikin forward and upward.

The jaw thrust maneuver is the safest approach to opening the airway of a victim who has a possible head, neck, or spinal injury because it can be accomplished by minimizing neck movement.



### Professional Adult Jaw Thrust, Single

PP-JTM-100M-MS  
PP-JTM-100M-DS



### Professional Adult Jaw Thrust, 4-pack

PP-JTM-100M-MS  
PP-JTM-100M-DS



- Multi-packs (pg 5) come with 10 Face-Shield/Lung-Bags per manikin and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin



# PROFESSIONAL

## PRESTAN Instructor Kit



The **Professional Instructor Kit** combines our most popular CPR training manikins (one Professional Adult, one Professional Female & two Professional Infants) along with two AED UltraTrainers® and a collection of the most commonly used CPR Training Supplies in one single collection, packaged in a heavy-duty canvas bag with wheels & telescoping handle for enhanced portability.



### Professional Instructor Kit

PP-INSTRKIT-1-MS  
PP-INSTRKIT-1-DS  
PP-INSTRKIT-1-MSDS



- Instructor Kit comes with (2) Adult CPR Manikin Training Shirts, (2) Infant CPR Manikin Training Shirts, (50) Adult Face-Shield/Lung-Bags, (50) Infant Face-Shield/Lung-Bags, (2) Adult Bag Valve Masks, (2) Infant Bag Valve Masks, (4) PRESTAN Kneeling Pads, (1) Carry Bag w/ Roller Wheels & Telescopic Handle



# PROFESSIONAL / PORTABLE

## PRESTAN Specialty Packs



The **PRESTAN Specialty Packs** combine our most popular products packaged together in single collections making them ideal as new instructor starter kits and for classroom or larger group training scenarios.

The **TAKE2® Kit** includes two Adult Manikins, two Infant Manikins, and two AED UltraTrainers—an ideal classroom starter kit for new instructors. Also available with one Adult and one Female Manikin, two Infant Manikins, and two AED UltraTrainers in medium or dark skin tones.



### TAKE2 Professional Kit

PP-FM-600M-MS  
PP-FM-600M-DS  
PP-FM-600M-MSDS



### TAKE2 Professional Kit

PP-AM-601M-MS  
PP-AM-601M-DS



### TAKE2 Ultralite Kit

PP-ULFM-600M-MS  
PP-ULFM-600M-DS  
PP-ULFM-600M-MSDS



- Multi-packs (pg 5) come with 10 Face-Shield/Lung-Bags per manikin and a nylon carrying case with a detachable shoulder strap
- Manikin shirts included with every manikin.



# PORTABLE

## PRESTAN Ultralite® Manikin



The **Ultralite Manikin** is our most compact, lightweight, portable manikin available in single, 4-pack and 12-pack configurations making it ideal for virtual and on-the-go training.

This durable manikin offers an affordable solution for CPR training, with all the quality and realism you expect from PRESTAN.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available with or without CPR feedback and in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty

The **Ultralite Piston with CPR Feedback** provides instructors with an option to upgrade their current non-monitored PRESTAN Ultralite Manikins with CPR rate and depth feedback.




### Ultralite Piston with CPR Feedback

RPP-ULPISTONM-1-MS (single)  
RPP-ULPISTONM-4-MS (4-pack)



### Ultralite Manikin, Single

PP-ULM-100M-MS  
PP-ULM-100M-DS



- Manikin shirts included with every manikin.

# PORTABLE

## PRESTAN Ultralite® Manikin



The Ultralite Manikin's unique, clamshell design with stackable torsos and heads creates a lightweight, compact solution ideal for transporting multiple units in packs of 4 (13 pounds) or 12 (39 pounds).

The **Ultralite Manikin, 4-pack** comes in an easy-to-carry bag with a detachable shoulder strap while the **Ultralite Manikin, 12-pack** offers a full classroom quantity in a portable roller bag with wheels and a telescoping handle for effortless maneuverability making it easy for lifting in and out of vehicles.



### Ultralite Manikin, 4-pack

PP-ULM-400M-MS  
PP-ULM-400M-DS



### Ultralite Manikin Diversity Kit, 4-pack

PP-ULM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)



### Ultralite Manikin, 12-pack

PP-ULM-1200M-MS  
PP-ULM-1200M-DS



### Ultralite Manikin Diversity Kit, 12-pack

PP-ULM-1200M-MSDS  
(6 medium skin tone & 6 dark skin tone)

- Manikin shirts included with every manikin.



# PORTABLE

## PRESTAN Infant Ultralite® Manikin



As with all PRESTAN manikins, the **Infant Ultralite Manikin** provides a realistic training experience with anatomically correct features, accurate head tilt and visible chest rise while providing intuitive feedback during CPR training.

This durable manikin offers a lightweight, affordable solution for CPR training, with all the quality and realism you expect from PRESTAN.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



**Infant Ultralite Manikin, Single**

PP-IULM-100M-MS  
PP-IULM-100M-DS



**Infant Ultralite Manikin, 4-pack**

PP-IULM-400M-MS  
PP-IULM-400M-DS



**Infant Ultralite Manikin Diversity Kit, 4-pack**

PP-IULM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)



**Infant Ultralite Manikin, 12-pack**

PP-IULM-1200M-MS  
PP-IULM-1200M-DS



**Infant Ultralite Manikin Diversity Kit, 12-pack**

PP-IULM-1200M-MSDS  
(6 medium skin tone & 6 dark skin tone)

- Manikin shirts included with every manikin.

# PORTABLE

## PRESTAN Infant Ultralite® Manikin



The compact, lightweight design of the Infant Ultralite Manikin allows for easy transport and cost-effective shipping making it an ideal solution for on-the-go and virtual training.



### Infant Ultralite Manikin, 4-pack

Lightweight and easy to transport in 5.8-pound carry bag.



### Infant Ultralite Manikin, Single

Just 1.6-pounds with easy to carry drawstring bag.



### Infant Ultralite Manikin, 12-pack

Full classroom quantity in one 15.5 pound portable bag.

- Manikin shirts included with every manikin.



# PROFESSIONAL

## PRESTAN Professional AED Trainer PLUS



The **Professional AED Trainer PLUS** is a leading-edge AED Trainer that most closely resembles brand AED Trainers in size and weight providing a realistic training experience.

This AED trainer features customizable options to match various training environments, including a child button.

- Clear, confident voice prompts with optional CPR prompts to help the student keep pace with chest compressions
- Optional remote available; a single remote controls several AED Trainers within range
- Ability to toggle between 15:2 and 30:2 compression-ventilation ratio
- Easy to replace module containing the most recent industry guidelines in two languages - comes standard with English/Spanish or English/French; other language options available upon request
- For use with patented PRESTAN AED Trainer Pad Sensing System
- 3-year Limited Warranty



### Professional AED Trainer PLUS

PP-AEDT2-101 (single; English/Spanish)

PP-AEDT2-102 (single; English/French)

PP-AEDT2-401 (4-pack; English/Spanish)

PP-AEDT2-402 (4-pack; English/French)



### Professional AED Trainer PLUS Kit

PP-AEDT2-KIT-101

(single with Remote; English/Spanish)



# PORTABLE

## PRESTAN AED UltraTrainer®



The **AED UltraTrainer** is our most compact, lightweight AED trainer solution with full customization options to accommodate any classroom training situation.

This AED trainer resembles a live training experience like no other on the market and offers an affordable solution with all the quality and durability you expect from PRESTAN.

- Instructors can customize their training classes with options to turn on or off the compression metronome and/or the “give breaths” voice prompt
- Programming can toggle between two languages and is easily upgradeable with the most recent industry guidelines - comes standard with English/Spanish or English/French; other language options available upon request
- Ability to toggle between 15:2 and 30:2 compression-ventilation ratio
- For use with patented PRESTAN AED Trainer Pad Sensing System
- 3-year Limited Warranty



### AED UltraTrainer

- PP-AEDUT-101 (single; English/Spanish)
- PP-AEDUT-102 (single; English/French)
- PP-AEDUT-401 (4-pack; English/Spanish)
- PP-AEDUT-402 (4-pack; English/French)






# PRESTAN AED Trainer Pad Sensing System

The patented AED Trainer Pad Sensing System automatically detects when each AED trainer pad is placed on the manikin. New and improved AED Trainer Pads feature enhanced sensor performance for more reliable pad detection during training. Pads can now be inserted in either orientation for faster, easier replacement.

These high-quality pads last for over 100 applications. The upgraded adhesive adheres well to any manikin and removes cleanly without leaving residue—no more damaged pads or cleanup between sessions. Pad replacement remains simple and affordable—just replace the pad, not the cord.

Programming runs through a sensor in the center of the pad, simulating a live AED for a more realistic training experience. No remote control is necessary, although one is available for the Professional AED Trainer PLUS if desired.



Universal compatibility with both the PRESTAN Professional AED Trainer PLUS and the PRESTAN AED UltraTrainer.



## Universal AED Trainer Pads

PP-UTPAD-1 (single; 2 pads total)

PP-UTPAD-4 (4-pack; 8 pads total)



# PRESTAN CPR Training Supplies



**PRESTAN Kneeling Pads** provide cushioned knee and back comfort for students during CPR training classes.




**Kneeling Pads**  
PP-KPAD-4 (4-pack)  
PP-KPAD-24 (24-pack)




The patent-pending **PRESTAN Female Accessory** has been designed to *build confidence* in CPR training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration.

The PRESTAN Female Accessory simply clips onto the manikin torso and securely fits both the PRESTAN Professional Adult/Series 2000 & Ultralite Manikins.



## PRESTAN Female Accessory

- PP-FA-1-DS (dark; single)
- PP-FA-4-DS (dark; 4-pack)
- PP-FA-1-MS (medium; single)
- PP-FA-4-MS (medium; 4-pack)





# PRESTAN CPR Training Supplies



The **PRESTAN CPR Training Face Masks** are made of a clear, non-collapsible medical grade PVC constructed with a moldable bladder for proper seal on a range of manikins; for use with training BVM's. Available in 10-packs; with or without individually wrapped **PRESTAN Training Face Mask Adaptors**.



## CPR Training Face Mask

PP-AMASK-10 (adult/child; 10-pack)  
PP-IMASK-10 (infant; 10-pack)

## CPR Training Face Mask & Adaptor Bundle

PP-AMASK-ADT-10 (adult/child; 10-pack)  
PP-IMASK-ADT-10 (infant; 10-pack)



## CPR Training Face Mask Adaptor

10076-PPA (10-pack; individually wrapped)  
10076-PPA-50 (50-pack; individually wrapped)



# PRESTAN CPR Training Manikin Consumables



The **PRESTAN Lung-Bag/Face-Shields** are easy to install/remove and allow for visual chest rise to ensure students are receiving the most accurate, realistic CPR training experience possible.

The **PRESTAN Professional Adult Series 2000 Manikin Lung-Bag/Face-Shields** use the same easy to install PRESTAN design with the addition of a breath sensor to help monitor ventilation feedback along with the PRESTAN CPR Feedback App.

**PRESTAN Face Shields** allow for multiple students training on the same manikin without the need to change lung bags and are specifically designed to snap on and securely fit the face of each PRESTAN manikin.



**Professional Adult/Female Manikin**  
PP-ALB-50 & PP-AFS-50



**PRO+ Adult Manikin**  
PP-APLB-50 & PP-APFS-50



**Professional Child Manikin**  
PP-CLB-50 & PP-CFS-50



**Professional Infant Manikin**  
PP-IFS-50 & PP-ILB-50



The patented PRESTAN manikin clamshell torso allows for quick and easy lung bag insertion and removal through the mouth and out the neck.



**PRO+ Series 2000 Manikin**  
PP-VAPLB-50



**Professional Adult/Female Series 2000 Manikin**  
PP-VALB-50 & PP-AFS-50



**Series 2000 Child Manikin**  
PP-VCLB-50



**Series 2000 Infant Manikin**  
PP-VILB-50



**Ultralite Manikin**  
PP-ULB-50 & PP-UFS-50



**Infant Ultralite Manikin**  
PP-IFS-50 & PP-ILB-50









[www.prestan.com](http://www.prestan.com)