



PRESTAN's Most Advanced CPR Training Technology

The PRESTAN Series 2000 Complete Offering Offering *Pre-Order Now!*



PRESTAN's Professional Series 2000 manikin offering is expanding — and we're bringing you the most advanced CPR feedback technology yet. Available mid-April in Adult, PRO+, Female, Child, and Infant models, the Series 2000 offering provides real-time, high-quality advanced feedback to enhance CPR training and build student confidence.



Whether teaching students to respond to an adult, a child, or an infant, the Series 2000 Complete Offering gives you everything you need to ensure they feel prepared to take immediate action. With intuitive feedback to help teach high-quality CPR, including compression rate, depth, recoil, and ventilation — your students will develop the skills and confidence to deliver high-quality CPR and critical life-saving care.



Key Features:

- **Compliant With Industry Guidelines:** Fully compliant with current industry guidelines, including the most recent AHA Integrated Feedback Directive. The Series 2000 is also fully compliant with the Heart and Stroke Foundation of Canada's February 2021 *Use of Feedback Devices in Resuscitation Training Courses Policy*.
- **Comprehensive Offering:** Advanced CPR Feedback available for Adult, Female, Child, and Infant manikins — in medium or dark skin tones and in single or 4-pack configurations.
- **Dual Purpose CPR Training Solution:** In combination with PRESTAN's traditional CPR feedback light system, the Series 2000 manikins provide an option to train with either advanced or basic feedback.
- **Advanced Feedback:** Compression rate, depth, recoil, ventilation, hands-off time, and chest compression fraction are all displayed in real-time via our redesigned PRESTAN CPR Feedback app.
- **Redesigned App:** The new PRESTAN CPR Feedback app, available on Apple and Android devices, provides an improved user experience with easier setup, intuitive CPR feedback — with no need to calibrate breaths, and an expanded, growing list of ~30 languages supported.
- **Instructor-Friendly Tracking:** Monitor CPR feedback for up to six students simultaneously, with detailed data available for performance review and reporting.
- **Made in the USA:** Unlike any other branded or privately labeled CPR training manikin available worldwide, all PRESTAN CPR training manikins are made in the United States and come with a three-year limited warranty.

Be Among the First to Train with the Series 2000 Complete Offering

Prepare your students with the confidence and skills required to perform high-quality CPR when needed.

The PRESTAN Series 2000 Complete Offering will be available in mid-April 2025. Contact your [PRESTAN Authorized Distributor](#) to place your pre-order today!

Thank you for your commitment to enhancing CPR training and empowering others to save lives!



PRESTAN Products, LLC | www.prestan.com



