

ENDLESS TRAINING POSSIBILITIES

PRESTAN Series 2000 Manikins

The **Series 2000 Manikins** pair realistic form with next-level digital feedback. Visual chest rise and comprehensive performance metrics—rate, depth, recoil, ventilation quality, hands-off time, and chest compression fraction—stream directly to the PRESTAN CPR Feedback App via Bluetooth®.



EVERYTHING YOU NEED TO TRAIN WITH CONFIDENCE

- Train a wide range of learners with versatile adult, female, child, and infant Series 2000 configurations.
- Build better CPR skills with realistic feedback that helps learners understand and improve their performance.
- Gain actionable training insights with intuitive feedback technology designed to support instructors and learners.
- Train with confidence using durable manikins built for repeated, hands-on practice.
- Take training anywhere with portable designs that make setup quick and easy across classrooms, training centers, workplaces, and more.



PRESTAN CPR Feedback App: Download Now on Apple or Android devices

