

RAIZOR TRAINER



COMPACT



LIGHTWEIGHT



PORTABLE

Designed for simple, cost-effective AED training, the PRESTAN AED Raizor Trainer delivers everything instructors need—without unnecessary complexity.

Its compact, all-blue design fits comfortably in the hand and supports straightforward instruction with one shock training scenario, a fixed volume level, and clear English voice prompts.





AED TRAINERS

PRESTAN offers a complete lineup of AED Trainers, giving instructors the flexibility to choose the right features, size, and functionality for any training environment.



	PROFESSIONAL AED TRAINER PLUS	AED ULTRATRAINER	RAIZOR TRAINER
Size: L x W x H (In.)	2.75 x 7.25 x 10	7.25 x 5.375 x 2	5 x 3.75 x 1.5
Hard Cover Lid	✓		
Optional Remote	✓		
Pad Sensors	✓	✓	
Child Button	✓	✓	
Carry Bag	✓	✓	
CPR Metronome	✓	✓	
Pad Storage Tray	✓	✓	
Pad Storage Envelope			✓
# of Training Scenarios	5	5	1
Languages on Single Unit	2	2	1
Battery Size	(4) C	(3) C	(2) AA
Warranty (Years)	3	3	1
MSRP	SSS	SS	S