



PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

How to translate PRESTAN Date Codes:

3-digit Date Code: WWY

- First two digits are week and last digit is year.
- For example, a 3-digit date code of 302 is 30th week of 2022 -or- the week of July 28th, 2022.

4-digit Date Code: YDDD

- First digit is year and the last three digits are the Julian date.
- For example, a 4-digit code of 2164 would be the 164th day of 2022 -or- June 13th, 2022.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Professional Adult / Female / Child Manikins:
Date code located on back of manikin torso.






PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Professional Adult / Female / Child Rate Monitor:
Date code located on inside of battery compartment.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Professional Infant Manikin:

Date code located on side of manikin torso under skin.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Ultralite Manikin:

Date code located on bottom lip/inside of manikin torso.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Ultralite Manikin Feedback Piston:

Date code located on inside of battery compartment.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Infant Ultralite Manikin:

Date code located on back/inside of diaper/manikin torso.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

AED UltraTrainer:

Date code located on inside of battery compartment.





PRESTAN®

PRESTAN Manufacturing Date Code Warranty Help Guide

Professional AED Trainer PLUS:

Date code located on inside of battery compartment.

