



# *Building Confidence* with PRESTAN CPR Manikins and AED Trainers



**PRODUCT  
CATALOG**

[www.prestan.com](http://www.prestan.com)



# Building Confidence with Feedback...





PRESTAN offers a full range of CPR Training Manikins and AED Trainers designed to meet the needs of all training environments - including **Portable**, **Professional**, and **Advanced** training scenarios. Whether an instructor is starting a new business or looking to upgrade their equipment, PRESTAN offers the most robust product offering in the market giving students the *confidence* and skills necessary to be ready for all lifesaving moments.



PRESTAN CPR Training Manikins have provided a feedback mechanism since their inception. All **PRESTAN Professional Manikins** and **PRESTAN Ultralite Manikins with CPR Feedback** come equipped with a CPR Monitor that provides both visual and audible feedback, helping students learn the proper compression rate and depth and allowing them to *build confidence* and skills needed to perform high quality CPR.

The **PRESTAN Professional Adult Series 2000 Manikin** provides advanced feedback through the addition of a Bluetooth connection to the **PRESTAN CPR Feedback App** which monitors rate, depth, recoil, ventilation, hands-off time, and chest compression fraction.



Both the **PRESTAN Professional AED Trainer PLUS** and the **PRESTAN AED UltraTrainer** provide a realistic training experience with clear, confident voice instructions and easy-to-follow CPR coaching feedback.





# Building Confidence with PRESTAN...



PRESTAN is committed to *building confidence* in CPR/AED training instructors and their students.

PRESTAN offers a full range of CPR Training Manikins and AED Trainers which are intuitive, durable, cost-effective, and provide immediate feedback.

PRESTAN also offers a full range of Training Supplies and Accessories including replacement lung bag/face shields, stand-alone face shields, training face masks & adaptors, manikin clothing, and kneeling pads.



## Find a PRESTAN Authorized Distributor

PRESTAN products are sold exclusively through a global network of PRESTAN Authorized Distributors.

Please visit the 'Find a Distributor' page on the PRESTAN website to locate a PRESTAN Authorized Distributor near you.



## Visit the PRESTAN Online Resource Center

PRESTAN provides an extensive set of educational tools and training resources that will give you and your students the *confidence* to be ready.

The PRESTAN Resource Center contains Product Instruction Sheets, Frequently Asked Questions, Troubleshooting & Product Help Guides, Instructional How-To & Product Overview Videos, and much more!





# Building Confidence with PRESTAN...



PRESTAN *builds confidence* in instructors by offering an industry-leading **3-year Limited Warranty** on all PRESTAN CPR Training Manikins and AED Trainers.

To learn more about the PRESTAN Limited Warranty, please contact a PRESTAN Authorized Distributor.



PRESTAN Products, LLC is headquartered in the United States and operates in the *heart* of the country near Cleveland, Ohio.

Unlike any other branded or private labelled CPR Training Manikin available throughout the world, all PRESTAN CPR Training Manikins are Made in the USA.

The company was founded in 2004 and has grown from a small start-up to a global company with warehousing facilities located throughout the world.

To learn more about PRESTAN Products, please visit [www.prestan.com](http://www.prestan.com).





# ADVANCED

## PRESTAN Professional Adult Series 2000



In combination with PRESTAN's traditional CPR feedback light system, the **Professional Adult Series 2000 Manikin** offers a unique, dual-purpose training solution providing instructors with an option to train with either basic or advanced CPR feedback including compression rate, depth, recoil, ventilation, hands-off time, and chest compression fraction via a Bluetooth enabled PRESTAN CPR Feedback App.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



\*Mobile device not included

The **Professional Adult Series 2000 Lung Bags** use the same easy to install PRESTAN design with the addition of a breath sensor. The Series 2000 lung bags are required for ventilation feedback via the PRESTAN CPR Feedback app.



**Series 2000 Lung Bags**  
PP-VALB-50 (50-pack)

### Professional Adult Series 2000, Single

PP-AM-2000-1-MS  
PP-AM-2000-1-DS



\*Mobile device not included


### Professional Adult Series 2000, 4-pack

PP-AM-2000-4-MS  
PP-AM-2000-4-DS

The **Professional Adult Series 2000 Manikin Add-On Kit** provides instructors with an option to upgrade their current PRESTAN Professional Adult and Jaw Thrust Manikins (*2018 models and newer*) with Advanced Feedback parameters.

### Series 2000 Add-On Kit

RPP-AAOK2000-1 (single)  
RPP-AAOK2000-4 (4-pack)



The **PRESTAN CPR Feedback App** provides high-quality, real-time, advanced CPR feedback via a Bluetooth connection with the PRESTAN Professional Adult Series 2000 Manikin.



**Download Now**  
on Apple or Android devices





# PROFESSIONAL PRESTAN Professional Adult Manikin



Realistic to the eye and the touch, the **Professional Adult Manikin** is unlike any other CPR training manikin on the market. Designed to provide each student with a real-life training experience along with immediate compression rate and depth feedback, PRESTAN manikins help instructors *build confidence* in their students when teaching high-quality CPR skills.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



## Professional Adult, Single

PP-AM-100M-MS

PP-AM-100M-DS



## Professional Adult, 4-pack

PP-AM-400M-MS

PP-AM-400M-DS



## Professional Adult Diversity Kit

PP-AM-400M-MSDS

(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap





# PROFESSIONAL PRESTAN Professional Female Manikin



The **Professional Female Manikin** is designed to transform CPR/AED training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration. Using the same trusted design as the Professional Adult Manikin, the Female Manikin is realistic to the eye and the touch, unlike any other CPR Training Manikin on the market.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



When used in combination with the PRESTAN Female Accessory (to represent a bra), Manikin Shirt, and AED Trainer, PRESTAN offers the most realistic training experience to help *build confidence* when learning these critical lifesaving skills.

## Professional Female, Single

PP-AFM-100M-MS

PP-AFM-100M-DS



## Professional Female, 4-pack

PP-AFM-400M-MS

PP-AFM-400M-DS



## Professional Female Diversity Kit

PP-AFM-400M-MSDS

(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap



# PROFESSIONAL PRESTAN PRO+ Manikin



The **PRO+ Adult Manikin** features a fully functioning oral and nasal airway along with an improved lung bag design for fast and easy installation. With the introduction of the PRO+ Adult Head design, PRESTAN is enhancing realism with a functional nasal airway and the ability to use a NARCAN® or NPA training device.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



When lubricated properly, a nasopharyngeal airway (NPA) can be inserted through the nose of the PRO+ Adult Head to simulate opening the nasal airway during training.



## PRO+ Adult, Single

PP-AMP-100M-MS  
PP-AMP-100M-DS

The **PRO+ Adult Head** is also sold separately, allowing you to upgrade any Professional Adult Manikin.



## PRO+ Adult Head

PP-APHEAD-1-MS  
PP-APHEAD-1-DS  
PP-APHEAD-4-MS  
PP-APHEAD-4-DS



## PRO+ Adult, 4-pack & Diversity Kit

PP-AMP-400M-MS  
PP-AMP-400M-DS  
PP-AMP-400M-MSDS



## PRO+ Lung Bag

PP-APLB-50 (50-pack)



## PRO+ Adult Face Shield

PP-APFS-50 (50-pack)



The **PRO+ Adult Face Shield** has been designed to fully cover both oral and nasal airways; the standard PRESTAN Professional Adult Face Shield can also be used.

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap





# PROFESSIONAL

## PRESTAN Professional Child/Infant Manikins



The **Professional Child Manikin** is designed to look like a child versus an adult with less prominent body structures and softer, more childlike facial features. It is smaller and thinner than the Professional Adult Manikin yet still realistic to the eye and touch.



### Professional Child, Single

PP-CM-100M-MS  
PP-CM-100M-DS



### Professional Child, 4-Pack

PP-CM-400M-MS  
PP-CM-400M-DS



### Professional Child Diversity Kit

PP-CM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)

The **Professional Infant Manikin** life-like features and weight allows for a realistic infant CPR training experience. Sculpted body contours and realistic skin feel help prepare students for real-life rescues with a head designed to simulate the way an actual infant's head moves.



### Professional Infant, Single

PP-IM-100M-MS  
PP-IM-100M-DS



### Professional Infant, 4-Pack

PP-IM-400M-MS  
PP-IM-400M-DS



### Professional Infant Diversity Kit

PP-IM-400M-MSDS  
(2 medium skin tone & 2 dark skin tone)

- Single packs come with 10 Face-Shield/Lung-Bags and a nylon carrying case
- 4-packs come with 50 Face-Shield/Lung-Bags and a nylon carrying case with a detachable shoulder strap



# PROFESSIONAL PRESTAN Specialty Packs



The **Professional Family Pack** combines two Professional Adult, one Child, and two Infant Manikins in one single collection.



## Professional Family Pack

PP-FM-500M-MS  
PP-FM-500M-DS

The **Professional Collection** combines a Professional Adult, Child, and Infant Manikin in one single collection.



## Professional Collection

PP-FM-300M-MS  
PP-FM-300M-DS


The **Professional Adult Jaw Thrust Manikin** gives instructors the option to teach with the jaw thrust or head-tilt/chin-lift maneuver allowing the student to grasp the manikin's lower jaw and lift with both hands displacing the jaw of the manikin forward and upward.

The jaw thrust maneuver is the safest approach to opening the airway of a victim who has a possible head, neck, or spinal injury because it can be accomplished by minimizing neck movement.



## Professional Adult Jaw Thrust, Single

PP-JTM-100M-MS  
PP-JTM-100M-DS



## Professional Adult Jaw Thrust, 4-pack

PP-JTM-100M-MS  
PP-JTM-100M-DS

- Multi-packs (pg 5) come with 10 Face-Shield/Lung-Bags per manikin and a nylon carrying case with a detachable shoulder strap





# PROFESSIONAL PRESTAN Instructor Kit



The **Professional Instructor Kit** combines our most popular CPR training manikins (two Professional Adults w/ Female Accessory & two Professional Infants) along with two AED UltraTrainers® and a collection of the most commonly used CPR Training Supplies in one single collection, packaged in a heavy-duty canvas bag with wheels & telescoping handle for enhanced portability.

## Professional Instructor Kit

PP-INSTRKIT-1-MS  
PP-INSTRKIT-1-DS  
PP-INSTRKIT-1-MSDS



- Instructor Kit comes with (2) Adult CPR Manikin Training Shirts, (2) Infant CPR Manikin Training Shirts, (50) Adult Face-Shield/Lung-Bags, (50) Infant Face-Shield/Lung-Bags, (2) Adult Bag Valve Masks, (2) Infant Bag Valve Masks, (4) PRESTAN Kneeling Pads, (1) Carry Bag w/ Roller Wheels & Telescopic Handle



# PROFESSIONAL / PORTABLE PRESTAN Specialty Packs



The **PRESTAN Specialty Packs** combine our most popular products packaged together in single collections making them ideal as new instructor starter kits and for classroom or larger group training scenarios.

The **TAKE2® Kits** combine two Adult Manikins, two Infant Manikins, and two AED UltraTrainers, making this collection an ideal classroom 'starter kit' for any new instructor.



### TAKE2 Professional Kit

PP-FM-600M-MS  
PP-FM-600M-DS  
PP-FM-600M-MSDS



### TAKE2 Ultralite Kit

PP-ULFM-600M-MS  
PP-ULFM-600M-DS  
PP-ULFM-600M-MSDS



- Multi-packs (pg 5) come with 10 Face-Shield/Lung-Bags per manikin and a nylon carrying case with a detachable shoulder strap





# PORTABLE PRESTAN Ultralite® Manikin



The **Ultralite Manikin** is our most compact, lightweight, portable manikin available in single, 4-pack and 12-pack configurations making it ideal for virtual and on-the-go training.

This durable manikin offers an affordable solution for CPR training, with all the quality and realism you expect from PRESTAN.

- Fully compliant with all worldwide curriculum & industry guidelines (when purchased or upgraded with CPR Feedback piston)
- Lightweight with fast & simple setup; easy to clean & maintain
- Available with or without CPR feedback and in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty

The **Ultralite Piston with CPR Feedback** provides instructors with an option to upgrade their current non-monitored PRESTAN Ultralite Manikins with CPR rate and depth feedback.





### Ultralite Manikin, Single

- PP-ULM-100M-MS (with feedback)
- PP-ULM-100M-DS (with feedback)
- PP-ULM-100-MS (w/out feedback)
- PP-ULM-100-DS (w/out feedback)

### Ultralite Piston with CPR Feedback

- RPP-ULPISTONM-1-MS (single)
- RPP-ULPISTONM-4-MS (4-pack)



# PORTABLE

## PRESTAN Ultralite® Manikin



The Ultralite Manikin's unique, clamshell design with stackable torsos and heads creates a lightweight, compact solution ideal for transporting multiple units in packs of 4 (13 pounds) or 12 (39 pounds).

The **Ultralite Manikin, 4-pack** comes in an easy-to-carry bag with a detachable shoulder strap while the **Ultralite Manikin, 12-pack** offers a full classroom quantity in a portable roller bag with wheels and a telescoping handle for effortless maneuverability making it easy for lifting in and out of vehicles.



### Ultralite Manikin, 4-pack

- PP-ULM-400M-MS (with feedback)
- PP-ULM-400M-DS (with feedback)
- PP-ULM-400-MS (w/out feedback)
- PP-ULM-400-DS (w/out feedback)



### Ultralite Manikin Diversity Kit, 4-pack

- PP-ULM-400M-MSDS (with feedback)
- PP-ULM-400-MSDS (w/out feedback)
- (2 medium skin tone & 2 dark skin tone)



### Ultralite Manikin, 12-pack

- PP-ULM-1200M-MS (with feedback)
- PP-ULM-1200M-DS (with feedback)
- PP-ULM-1200-MS (w/out feedback)
- PP-ULM-1200-DS (w/out feedback)



### Ultralite Manikin Diversity Kit, 12-pack

- PP-ULM-1200M-MSDS (with feedback)
- PP-ULM-1200-MSDS (w/out feedback)
- (6 medium skin tone & 6 dark skin tone)





# PORTABLE

## PRESTAN Infant Ultralite® Manikin



As with all PRESTAN manikins, the **Infant Ultralite Manikin** provides a realistic training experience with anatomically correct features, accurate head tilt and visible chest rise while providing intuitive feedback during CPR training.

This durable manikin offers a lightweight, affordable solution for CPR training, with all the quality and realism you expect from PRESTAN.

- Fully compliant with all worldwide curriculum & industry guidelines
- Lightweight with fast & simple setup; easy to clean & maintain
- Available in both medium skin (MS) or dark skin (DS) tones
- All PRESTAN CPR Training Manikins are Made in the USA
- 3-year Limited Warranty



**Infant Ultralite Manikin, 4-pack**

PP-IULM-400M-MS

PP-IULM-400M-DS



**Infant Ultralite Manikin, Single**

PP-IULM-100M-MS

PP-IULM-100M-DS



**Infant Ultralite Manikin Diversity Kit, 4-pack**

PP-IULM-400M-MSDS

(2 medium skin tone & 2 dark skin tone)



**Infant Ultralite Manikin, 12-pack**

PP-IULM-1200M-MS

PP-IULM-1200M-DS



**Infant Ultralite Manikin Diversity Kit, 12-pack**

PP-IULM-1200M-MSDS

(6 medium skin tone & 6 dark skin tone)



# PORTABLE PRESTAN Infant Ultralite® Manikin



The compact, lightweight design of the Infant Ultralite Manikin allows for easy transport and cost-effective shipping making it an ideal solution for on-the-go and virtual training.



**Infant Ultralite Manikin, 4-pack**  
Lightweight and easy to transport  
in 5.8-pound carry bag.



**Infant Ultralite Manikin, Single**  
Just 1.6-pounds with easy to  
carry drawstring bag.



**Infant Ultralite Manikin, 12-pack**  
Full classroom quantity in one  
15.5 pound portable bag.





# PROFESSIONAL

## PRESTAN Professional AED Trainer PLUS



The **Professional AED Trainer PLUS** is a leading-edge AED Trainer that most closely resembles brand AED Trainers in size and weight providing a realistic training experience.

This AED trainer features customizable options to match various training environments, including a child button.

- Clear, confident voice prompts with optional CPR prompts to help the student keep pace with chest compressions
- Optional remote available; a single remote controls several AED Trainers within range
- Ability to toggle between 15:2 and 30:2 compression-ventilation ratio
- Easy to replace module containing the most recent industry guidelines in two languages - comes standard with English/Spanish or English/French; other language options available upon request
- For use with patented PRESTAN AED Trainer Pad Sensing System
- 3-year Limited Warranty



### Professional AED Trainer PLUS

- PP-AEDT2-101 (single; English/Spanish)
- PP-AEDT2-102 (single; English/French)
- PP-AEDT2-401 (4-pack; English/Spanish)
- PP-AEDT2-402 (4-pack; English/French)



### Professional AED Trainer PLUS Kit

- PP-AEDT2-KIT-101 (single with Remote; English/Spanish)



# PORTABLE PRESTAN AED UltraTrainer®



The **AED UltraTrainer** is our most compact, lightweight AED trainer solution with full customization options to accommodate any classroom training situation.

This AED trainer resembles a live training experience like no other on the market and offers an affordable solution with all the quality and durability you expect from PRESTAN.

- Instructors can customize their training classes with options to turn on or off the compression metronome and/or the “give breaths” voice prompt
- Programming can toggle between two languages and is easily upgradeable with the most recent industry guidelines - comes standard with English/Spanish or English/French; other language options available upon request
- Ability to toggle between 15:2 and 30:2 compression-ventilation ratio
- For use with patented PRESTAN AED Trainer Pad Sensing System
- 3-year Limited Warranty



## AED UltraTrainer

- PP-AEDUT-101 (single; English/Spanish)
- PP-AEDUT-102 (single; English/French)
- PP-AEDUT-401 (4-pack; English/Spanish)
- PP-AEDUT-402 (4-pack; English/French)






# PRESTAN AED Trainer Pad Sensing System



The patented **AED Trainer Pad Sensing System** automatically detects when each AED trainer pad is placed on the manikin.

These high-quality AED trainer pads last for over 100 applications! The pad adhesive removes from and adheres well to any manikin on the market without leaving messy residue behind - no more damaged pads or cleaning adhesive off your manikins between sessions. Pad replacement is affordable and easy - just replace the pad, not the cord.



Programming runs through a sensor in the center of the pad's adhesive which simulates a live AED, allowing for a more realistic training experience. No remote control is necessary, although one is available for the Professional AED Trainer PLUS if desired.



## Universal AED Trainer Pads

PP-UTPAD-1 (single; 2 pads total)

PP-UTPAD-4 (4-pack; 8 pads total)



# PRESTAN CPR Training Supplies



**PRESTAN CPR Training Shirts** provide instructors with the ability to *build confidence* during training with more real-life scenarios for their students.

The Adult/Child/Ultralite shirt (one size fits all) includes an easy-to-remove front zipper with a non-skid back. The Infant shirt features an easy-to-remove 'hook and loop' design.



**CPR Training Shirt  
Adult/Child/Ultralite Manikin**  
PP-ACSHIRT-4 (4-pack)



**CPR Training Shirt  
Infant Manikin**  
PP-ISHIRT-4 (4-pack)

**PRESTAN Kneeling Pads** provide cushioned knee and back comfort for students during CPR training classes.



**Kneeling Pads**  
PP-KPAD-4 (4-pack)  
PP-KPAD-24 (24-pack)






# PRESTAN CPR Training Supplies



The patent-pending **PRESTAN Female Accessory™** has been designed to *build confidence* in CPR training, improve CPR survival rates in women, encourage dialogue & breakdown barriers, and eliminate gender bias in CPR administration.

The PRESTAN Female Accessory simply clips onto the manikin torso and securely fits both the PRESTAN Professional Adult/Series 2000 & Ultralite Manikins.



## PRESTAN Female Accessory

- PP-FA-1-DS (dark; single)
- PP-FA-4-DS (dark; 4-pack)
- PP-FA-1-MS (medium; single)
- PP-FA-4-MS (medium; 4-pack)



The **PRESTAN CPR Training Face Masks** are made of a clear, non-collapsible medical grade PVC constructed with a moldable bladder for proper seal on a range of manikins; for use with training BVM's. Available in 10-packs; with or without individually wrapped **PRESTAN Training Face Mask Adaptors**.



## CPR Training Face Mask

- PP-AMASK-10 (adult/child; 10-pack)
- PP-IMASK-10 (infant; 10-pack)



## CPR Training Face Mask & Adaptor Bundle

- PP-AMASK-ADT-10 (adult/child; 10-pack)
- PP-IMASK-ADT-10 (infant; 10-pack)



## CPR Training Face Mask Adaptor

- 10076-PPA (10-pack; individually wrapped)
- 10076-PPA-50 (50-pack; individually wrapped)



# PRESTAN CPR Training Manikin Consumables



The **PRESTAN Lung-Bag/Face-Shields** are easy to install/remove and allow for visual chest rise to ensure students are receiving the most accurate, realistic CPR training experience possible.

The **PRESTAN Professional Adult Series 2000 Manikin Lung-Bag/Face-Shields** use the same easy to install PRESTAN design with the addition of a breath sensor to help monitor ventilation feedback along with the PRESTAN CPR Feedback App.

**PRESTAN Face Shields** allow for multiple students training on the same manikin without the need to change lung bags and are specifically designed to snap on and securely fit the face of each PRESTAN manikin.



**Professional Adult/Female Manikin**  
PP-ALB-50 & PP-AFS-50



**Professional Adult Series 2000 Manikin**  
PP-VALB-50 & PP-AFS-50



**PRO+ Adult Manikin**  
PP-APLB-50 & PP-APFS-50



**Professional Child Manikin**  
PP-CLB-50 & PP-CFS-50



**Professional Infant Manikin**  
PP-IFS-50 & PP-ILB-50



**Ultralite Manikin**  
PP-ULB-50 & PP-UFS-50



The patented PRESTAN manikin clamshell torso allows for quick and easy lung bag insertion and removal through the mouth and out the neck.



**Infant Ultralite Manikin**  
PP-IFS-50 & PP-ILB-50







*Building Confidence to save lives*



@prestanproducts



[www.prestan.com](http://www.prestan.com)