



PRESTAN®

Prestan AED Trainers

In the market for a new AED Trainer? Which one best suits your training needs?



Prestan Professional AED Trainer PLUS

Features

Prestan AED UltraTrainer

Similar to Full-Size AED Unit is 9.5" x 7.25" x 2.75" and weighs 2.2 lbs.	Size and Weight	Compact and Portable Unit is 7.25" x 5.25" x 2" and weighs 1.2 lbs.
Yes	CPR Feedback	Yes
Yes	Pad Sensing System (PSS)	Yes
Yes	Pre-connected Pads	Yes
5 Scenarios Available	Shock Simulation	5 Scenarios Available
Yes	Auto / Semi-Automatic Shock	Yes
Toggle On/Off	Rescue Breaths	Toggle On/Off
Toggle On/Off	Compression Metronome	Toggle On/Off
Yes	Child Button	Yes
Dual Adult/Child Pads	Pediatric Pads	Dual Adult/Child Pads
Available	Remote Control	No
Yes	Low Battery Prompt	Yes
3 Years	Prestan Warranty	3 Years

Visit www.prestanproducts.com for full details on these quality products.