



## All Prestan Professional Monitored Manikins are now equipped with an upgraded **BLUE** CPR Rate Monitor.

The BLUE Monitor provides the same great feedback, but now warns if the 120 compression upper limit is exceeded (as recommended in the 2015 AHA Guidelines).



**New Feature:** The flashing yellow light signals **SLOW DOWN!** When compressions exceed 120 per minute.



Monitor installs in diaper area of Infant Manikin (Monitor Part # RPP-IMON-1)

Monitor installs in back shoulder of Adult Manikin (Monitor Part # RPP-AMON-1) and Child Manikin (Monitor Part # RPP-IMON-1)



**PRESTAN**  
www.prestanproducts.com

### How it works:

When compressions are performed and reach the correct depth of 2"-2.4", an audible click will be heard to confirm.


A signal is then sent from the clicker to the monitor, with each compression, to instantly and continuously respond to compression speed with color-coded LED lighting... **Invaluable feedback for students and instructors!**

Rate and Depth Feedback is based on the recommendations in the 2015 AHA Guidelines. The original beige monitors are in compliance with the 2015 Guidelines without the 120 feature, but are no longer available from Prestan.


Blue Monitors may be purchased separately to easily upgrade your existing Prestan Professional Manikins, even if they did not have a monitor when purchased. **Contact your Prestan Authorized Distributor to order today!**

## All Prestan Professional Monitored Manikins are now equipped with an upgraded **BLUE** CPR Rate Monitor.

The BLUE Monitor provides the same great feedback, but now warns if the 120 compression upper limit is exceeded (as recommended in the 2015 AHA Guidelines).



**New Feature:** The flashing yellow light signals **SLOW DOWN!** When compressions exceed 120 per minute.



Monitor installs in diaper area of Infant Manikin (Monitor Part # RPP-IMON-1)

Monitor installs in back shoulder of Adult Manikin (Monitor Part # RPP-AMON-1) and Child Manikin (Monitor Part # RPP-IMON-1)



**PRESTAN**  
www.prestanproducts.com

### How it works:

When compressions are performed and reach the correct depth of 2"-2.4", an audible click will be heard to confirm.

A signal is then sent from the clicker to the monitor, with each compression, to instantly and continuously respond to compression speed with color-coded LED lighting... **Invaluable feedback for students and instructors!**

Rate and Depth Feedback is based on the recommendations in the 2015 AHA Guidelines. The original beige monitors are in compliance with the 2015 Guidelines without the 120 feature, but are no longer available from Prestan.

Blue Monitors may be purchased separately to easily upgrade your existing Prestan Professional Manikins, even if they did not have a monitor when purchased. **Contact your Prestan Authorized Distributor to order today!**