




PRESTAN®


PRESTAN AED UltraTrainer™

The **Prestan AED UltraTrainer™** provides the options instructors have been waiting for. Combining the most recent technology with realistic features not found on other AED Trainers, this product will excel in any training atmosphere. It is a compact and customizable unit providing quality, durability and affordability. The **Prestan AED UltraTrainer™** is sure to be an industry favorite.



Amazing Training Pads

- Pads are pre-connected and ready to use, unlike other AED trainers that use outdated plug-in technology
- Pads can be used 100 times (or more!) during CPR/AED classes
- One set of dual-graphics pads work for both adult and child training
- Pads have an embedded 'switch' which senses proper pad attachment to the manikin (Pad Sensing System) (patent pending)
- The silicone base adhesive on the pads, coupled with the foam structure, ensures easy application & removal



Units include one set of Adult/Child Pads, pre-installed batteries, instructions, extra connector buttons and a carrying case

Product Part Numbers

- ITEM # PP-AEDUT-101 Single AED UltraTrainer - English/Spanish
- ITEM # PP-AEDUT-401 Four AED UltraTrainers - English/Spanish
- ITEM # PP-AEDUT-102 Single AED UltraTrainer - English/French
- ITEM # PP-AEDUT-402 Four AED UltraTrainers - English/French

Features	Benefits	Prestan AED UltraTrainer™	Other Universal AED Trainers
AED Trainer Pads			
Pad Sensing System (PSS) (patent pending)	Senses when pads are placed on manikin then advances to analyzation which best simulates live AEDs. Others must plug in the pads cord to analyze which is technology no longer used with AEDs.	✓	✗
Pads Cord Preconnected	Simulates most modern live AEDs. Others must plug in the pads cord to analyze - technology no longer used with AEDs.	✓	✗
100+ Pad Application Life	Quality Pads last longer saving your budget. Other use low quality pads that need to be replaced more often.	✓	✗
CPR Coaching			
CPR Coaching, Ventilation (Give two breaths) Command	Giving breaths is part of CPR certification curriculums. Other AED trainers do not coach ventilation as required in CPR guidelines.	✓	✗
CPR Coaching, Compression Metronome	CPR Metronome beeps at pace recommend by CPR guidelines. Helps students keep proper compression rate.	✓	✓
On/Off Compression Metronome	Customize training class by turning off metronome if needed.	✓	✗
On/Off Ventilation Command	Customize training class by turning off ventilations if needed. Valuable when teaching Compression-Only CPR.	✓	✗
Child AED Training			
Child Button	Most accurately simulates modern AEDs resulting in a realistic training class.	✓	✗
Child Commands	Most accurately simulates modern AEDs resulting in a realistic training class.	✓	✗
Duel Graphic Pad	Simulates modern live AEDs using adult and child graphics on same pad. Others have separate Child pads adding expense.	✓	✗
Battery Life			
Long Life C Batteries	Class interruptions are minimized with long lasting C batteries which last 3x longer than AA batteries. Others use shorter life AA batteries.	✓	✗
Low Battery Warning	Notification helps avoid interruption of dead batteries during training class.	✓	✗
Shock Simulation			
Automatic or Semi-Automatic Shock Simulation	Customize AED training to the customer's AED.	✓	✗
Warranty			
3-Year Warranty	Durable product backed by Prestan. Others offer 1 year or less.	✓	✗