


## Introducing the Prestan Professional Infant Manikin – guaranteed to be the future of CPR Infant training manikins!




**PRESTAN**



- Life-like features and weight allows for realistic infant CPR training
- Sculpted body contours and realistic skin feel help prepare students for real life rescues.
- A breakthrough in design, Prestan Professional manikins let instructors set up and pack up at training sites in minutes, and provides students with a realistic training experience that prepares them for real rescues.
- Available individually or in convenient multipacks of four.
- A unique clamshell design that accommodates an easy-to-insert face shield lung bag

- Realistic-feeling skin that won't shrink or stretch, even with repeated cleanings
- A clicker mechanism that helps students learn to use the right amount of force in order to compress the victim's chest to the correct depth – great for boosting students' confidence in their ability to administer adequate compressions.
- A head that simulates the way an actual victim's head would move. Face/head tilt enhances the CPR learning experience and teaches students how to open the airway during CPR.
- Full chest rise and lung bag ensures that students receive the most accurate and realistic training.



- An optional CPR Rate Monitor that provides instant visual feedback on compression rates
  - One green light shows students they have reached adequate compression rates – two green lights tell students to maintain their current rate.
  - Realistic looking and acting manikin is more acceptable to both instructors and students.
  - Gives each student a realistic experience that provides for the best "real life" CPR outcomes.
- Easy-to-insert face shield lung bag speeds setup of the manikin.
  - Face shield is secured so it won't slip and there's no need for adhesives.

701 Beta Drive Suite 2 Mayfield Village, OH 44143  
(P) 440-229-5100 (F) 440-229-5111